



SURGICAL SOLUTIONS

WARMING  
SOLUTIONS



DESIGNED BY NURSES  
FOR NURSES™



- INNOVATIVE DESIGNS
- EFFORTLESS OPERATION
- AFFORDABLE SOLUTIONS

# WEG Warmer™ WARMING SYSTEMS

## CONTROL PANEL

- Easy to use
- Easy to read
- Turns on at 104°F every time  
(Following AORN Guidelines)

## BASIN

- Precise temperature control in the basin
- Patented Air Gap Detection System to prevent melt damage to the drape
- Rounded corners to protect the drape

## BASE

- Large casters ensure it's easy to move
- Low center of gravity to increase stability



## DRAPE

- Embedded frame ensures precise placement every time - helping to eliminate air gaps
- Dual layer drape-ensures sterility



## WEG WARMER™ WARMING DRAPES

Product Code	Description	Dimensions	Case Quantity
SD-REC44	Solution Warming Drape-Rectangular Frame	44"x44"	24
SD-REC60	Solution Warming Drape-Rectangular Frame	60"x44"	24
SD-CIR44	Solution Warming Drape-Circular Frame	44"x44"	24
SD-DISC48	Solution Warming Drape-Slush Disc	48"x48"	24
SD-DISC66	Solution Warming and Slush Drape Rectangular Frame-Slush Disc	66"x44"	24
SD-DISC66W	Solution Warming and Slush Drape Rectangular Frame-Slush Disc	66"x48"	24

## WEG WARMER™ WARMING DEVICES

Equipment Model	Description	Dimensions
SW-R38L	WEG Warmer™ Device 3.8 Liter Rectangular	11.5" x 5.25" x 4"
SW-R58L	WEG Warmer™ Device 5.8 Liter Rectangular with retractable scope cradle	19.5" x 5.0" x 4"
SW-C72L	WEG Warmer™ Device 7.2 Liter Circular	14.5" O.D. Round x 4"

Contact WEG Surgical Solutions at: [customerservice@wegsurgical.com](mailto:customerservice@wegsurgical.com)  
(888) 992-9888

For More Info: [www.wegsurgicalsolutions.com](http://www.wegsurgicalsolutions.com)